


PROVA PRATICA 1

Il candidato formuli ipotesi diagnostiche dalla visione dell’immagine


Overview Report SPECTRALIS® Tracking Laser Tomography

IR 30° ART + OCT 15° (8.3 mm) ART (9) Q: 24 [HR]




200 µm

IR 30° ART + OCT 15° (8.3 mm) ART (8) Q: 28 [HR]



200 µm

IR 30° ART + OCT 15° (8.3 mm) ART (9) Q: 29 [HR]



200 µm

Notes:

200 µm

200 µm

200 µm

Date:

Signature: